


體驗新一代核酸萃取的時機到了



任選買2送1
(7月31日前)

Spin DNA / RNA / Plasmid 純化套組

高純度

新結構管柱大幅
降低buffer及鹽
類殘留率

高濃縮

新式模層可超低
體積回溶，濃度
大幅提升

高產量

絕佳回收率，
頂尖產量看得見

低裂解

優化的步驟使核
酸完整性優於同
類產品

省荷包

Column及
Buffer等可單獨
購買補充，避免
浪費

愛地球

包裝堅持環境友
善材質，塑料用
量顯著低於市場
同類產品

Monarch[®] Spin RNA Isolation Kit (mini)

- 從各類生物樣本得到NGS等級RNA
- 適用肌肉, 腦, 植物等難處理樣本
- Elution 體積可低至10µl; 吸附上限100ug
- 也能用於回收短片段RNA(<200nt)
- 附gDNA column與DNase I, 清除gDNA污染
- 附樣本安定劑, 可在室溫保護樣本RNA

Monarch[®] Spin gDNA Extraction Kit

- 從各類生物樣本得到NGS等級DNA
- 大鼠組織測試結果產量明顯優於同類產品
- 多數樣本可獲得peak size >50kb之gDNA
- 附RNase A, RNA殘留極低(<0.1%)
- 如需Megabase-size DNA可選擇Monarch HMW DNA Extraction kit (T3050 / T3060)

Monarch[®] Spin DNA Gel Extraction Kit Spin PCR & DNA Cleanup Kit

- 新式管柱避免抑制物殘留，純度絕佳
- Elution體積可低至5µl, 濃度問題迎刃而解
- 不需調整pH值或添加isopropanol
- 縮短流程, cleanup僅需5分鐘
- 也能純化短片段DNA或oligonucleotide

Monarch[®] Spin Plasmid Miniprep Kit

- 快速從菌液中純化Plasmid DNA
- 以多種質體測試，濃度/產量/純度皆優於知名同類產品
- 顯色buffer系統，操作直覺便利不出錯

Order information:

Monarch Spin RNA Isolation Kit (Mini)	T2110S, 50 preps
Monarch Spin RNA Cleanup Kit (10 µg)	T2030S/L, 10/100 preps
Monarch Spin RNA Cleanup Kit (50 µg)	T2040S/L, 10/100 preps
Monarch Spin RNA Cleanup Kit (500 µg)	T2050S/L, 10/100 preps
Monarch Spin gDNA Extraction Kit	T3010S/L, 50/150 preps
Monarch HMW DNA Extraction Kit for Cell & Blood	T3050S/L, 5/50 preps

Monarch HMW DNA Extraction Kit for Tissue	T3060S/L, 5/50 preps
Monarch Spin Plasmid Miniprep Kit	T1110S/L, 50/250 preps
Monarch Spin DNA Gel Extraction Kit	T1120S/L, 50/250 preps
Monarch Spin PCR & DNA Cleanup Kit	T1130S/L, 50/250 preps
Monarch Mag Viral DNA/RNA Extraction Kit	T4010S/L/X, 100/600/1800 preps

